

$^1\text{H}\{^{31}\text{P}\}$ – Phosphorus Decoupled Proton NMR on the UI400 or U500

Decoupling ^{31}P from ^1H on the UI400 or U500 is easily accomplished with the standard QUAD probe. This handout outlines a procedure to acquire $^1\text{H}\{^{31}\text{P}\}$ spectra, where the phosphorus decoupling can be selective or non-selective. If you want to run this experiment on the U400 or VXR500 a cabling change is required. Check with Vera Mainz before attempting to run this experiment on any instruments EXCEPT the UI400 or U500.

Preliminaries:

- **Insert, Lock and Shim** as usual.
- If you are doing variable temperature work AND want to acquire $^1\text{H}\{^{31}\text{P}\}$ spectra, check with Lab staff as the temperature change might affect the efficiency of the decoupling.

Acquire Preliminary Spectra:

NOTE: Even if you've acquired this data before, you need to run them again.

jexp1 join experiment 1

Setup and run a standard ^1H experiment.

jexp2 join experiment 2

Setup and run a standard ^{31}P experiment. Make sure you can see all the ^{31}P peaks you plan to decouple.

1) If you want to decouple all ^{31}P peaks at the same time, place the left cursor at the center of the region of interest and type movetof <rt>. Record tof=_____.

2) If you plan to acquire selectively decoupled $^1\text{H}\{^{31}\text{P}\}$ spectra, place the cursor on the center of a ^{31}P peak of interest, type movetof. Record tof=_____.

If you have more than one ^{31}P peak to decouple, you need to reload the *FILE-31P-REFN* fid and wft for each frequency. After the wft, place the cursor on the center of the ^{31}P peak of interest, type movetof. Record tof=_____.

Continue to repeat this step for each ^{31}P peak of interest.

NOTE: You need to reload the reference ^{31}P spectrum each time before you execute another movetof command!!

$^1\text{H}\{^{31}\text{P}\}$ Decoupling Setup – decoupling all ^{31}P peaks at the same time:

jexp3 join experiment 1
mp(1,3) move the ^1H parameters from exp1 to exp3
xdec **sets up for $^1\text{H}\{^{31}\text{P}\}$ decoupling**
 Enter nucleus to be decoupled (e.g. P31): P31<rt>

Check that dn=P31

dof=_____

set the tof value from the ^{31}P spectrum,
above; if you have multiple values, you can
do them one at a time or in an array.

nt=1, if possible, but the minimum number of scans needed
ga

NOTE: To acquire a P31-coupled and decoupled array (so you can compare/plot them easily) set

dm='nnn', 'nny'

gain='y'

Turns off the autoscale gain function for
arrays

dssa or dssh

display arrayed data stacked vertically or horizontally

pl('all') pscale page

plot stacked data

$^1\text{H}\{^{31}\text{P}\}$ Decoupling Setup – decoupling ^{31}P peaks one at a time:

Use the same setup as above, with the following changes:

jexp3

join experiment 1

mp(1,3)

move the 1H parameters from exp1 to exp3

xdec

sets up for $^1\text{H}\{^{31}\text{P}\}$ decoupling

Enter nucleus to be decoupled (e.g. P31): P31<rt>

Check that dn=P31

dof=_____

set the tof value from the ^{31}P spectrum,
above; if you have multiple values, you can
do them one at a time or in an array.

dm='nny'

dmm='c'

dmf=250

dpwr=20 (This may change depending on how close your peaks are to each other and the JPP coupling constant. Check with Lab staff if you have problems.)

nt=1, if possible, but the minimum number of scans needed

ga

Example: Bis(diphenylphosphino)methane monoxide (dppmO) – $\text{Ph}_2\text{P}-\text{CH}_2-(\text{Ph}_2)\text{P}=\text{O}$ (data from *Inorg. Chem.* *v19(7)*, **1980**, 1982-1987):

$$^2J_{\text{P-CH}} = <0.5 \text{ Hz}; ^2J_{\text{P(O)CH}} = 12.5 \text{ Hz}$$


δ ^1H : $\text{CH}_2 = 3.06$ (doublet) (in my data, the doublet at $\sim \delta$ 3.22 is the CH_2 from dppmO; the triplet at $\sim \delta$ 3.74 is the CH_2 from dppmO₂)

δ ^{31}P : $\text{P}-\text{CH}_2 = -28.4$ (in my data, found at -26.7ppm)


δ ^{31}P : $\text{P}(\text{O}) \text{CH}_2 = 27.7$ (in my data, found at 33.3ppm)

δ ^{31}P dppmO₂: $\text{P}(\text{O}) \text{CH}_2 = 24.2$ (in my data, found at 28.6ppm)

dppmO – Ph₂P-CH₂-(Ph₂)P=O
All ³¹P Peaks Decoupled



³¹P{¹H} Spectrum dppmO and dppmO₂



dppmO – Ph₂P-CH₂-(Ph₂)P=O
¹H Spectra showing ³¹P Peaks Decoupled One at a Time

