

UVU422

MFL-15AUG91CD  
VVM-17APR92UD

### External Referencing of Samples on the U400

There are several ways to reference the chemical shifts of a given sample.

1. Reference to the solvent. This method is most commonly used for  $^1\text{H}$  and  $^{13}\text{C}$  NMR where an organic solvent such as  $\text{CDCl}_3$  is used which has residual protons and carbon-13 at natural abundance.
2. Reference to an added standard. Tetramethylsilane (TMS) is the standard reference for  $^1\text{H}$ ,  $^{13}\text{C}$  and  $^{29}\text{Si}$  NMR.
3. Reference to a compound in a capillary tube which is placed in the research sample. In this case, the standard might be insoluble in the lock solvent chosen or it might be reactive with the solvent or research sample.
4. Reference to an "external" standard. This usually means that a standard sample is inserted in the spectrometer, a frequency measurement taken, the standard removed, the research sample inserted, a spectrum obtained and the frequency calibration applied to the spectrum.

This fourth method of referencing is the focus of this handout. Three possible cases exist and are given in order of preference for accuracy of the reference.

1. The standard sample and the research sample are both in the same lock solvent at similar concentration and ionic strength
2. The standard sample is in a different lock solvent than the research sample.
3. The standard sample is not in a lock solvent, or the research sample is not in a lock solvent.

Three examples are given below.

**Example One:** The standard sample has a lock solvent.

$^{19}\text{F}$  research sample in  $\text{CDCl}_3$  and  $\text{CF}_3\text{C}_6\text{H}_5$  standard in  $\text{C}_6\text{D}_6$ . (Note that  $\text{CDCl}_3$  at 7.26 ppm and  $\text{C}_6\text{D}_6$  at 7.24 ppm are close enough in lock frequency that this will work.)

insert the research sample

```
jexp1 <>
LC Main Menu
LC Setup
LC Nucleus,Solvent
LC F19
LC CDCl3
rts('current shim library') <>
su <>
```

*use the standard setup for F19*

```
join experiment 1
activate main menu display
activate setup menu
select nucleus, solvent menu
select nucleus
select solvent
retrieve shims
perform experiment set up
```

When set up is complete:

load='n' <>  
dg <>

set load shim values to no  
display acquisition parameters

lock the sample and crudely shim Z1 and Z2 and disconnect the VNMR Acquisition window  
set nt to some reasonable number that will allow the peak(s) to be visible

ga <>

start acquisition (will wft when complete)

when acquisition is complete:

f full <>  
aph <>

display full sweep width to screen  
autophase

make sure that the peak(s) are visible at this sweep width and number of transients

LC Connect (in the VNMR acquisition window)

IMPORTANT:

LC Lock

If you do not turn the lock off before  
ejecting the sample, your reference

LC Lock off

will not be correct.

eject the research sample

insert the standard sample

*do not change any parameters*

lock on the standard by

LC Lock off

turn lock off

LC Lock on

turn lock on

*do not change Z0*

it might be necessary to increase lockpower and/or lockgain to lock on the standard

once locked, crudely shim Z1 and Z2 and disconnect the VNMR Acquisition Window

nt=4 <>

set number of transients

ga <>

start acquisition (will wft when complete)

when acquisition is complete:

f full <>  
aph <>

display full sweep width to screen  
autophase

LC HOLD, then release

move first cursor to the left of reference line

RC HOLD, then release

move second cursor to the right of reference line

LC Expand

expand region inside cursors

LC HOLD, then release

move first cursor to center of line

nl <>

select nearest line

rl(-63.73p) <>

set reference (in ppm)

f full <>

display full sweep width to screen

LC HOLD, then release

move first cursor far right of spectrum, to the last  
data point on the screen

LC Mark

request frequency

cr = \_\_\_\_\_

write down the frequency (which is in hertz)

**It is a good idea at this time to save the reference FID.**

LC Connect (in the VNMR acquisition window)

LC Lock

LC Lock off

eject the standard sample

insert the research sample

lock on the research sample without changing Z0

the research sample should now be shimmed as normal and the VNMR Acquisition window disconnected

*do not change any other parameters*

set nt so that the peak(s) will be visible

ga <>

IMPORTANT:

If you do not turn the lock off before

ejecting the sample, your reference

will not be correct.

start acquisition (will wft when complete)

When acquisition is complete:

f full <>

aph <>

LC HOLD, then release

rl(cr value in hertz) <>

f full <>

display full sweep width to screen

autophase

move first cursor far right of spectrum, to the last

data point on the screen

set reference (in hertz)

display full sweep width to screen

**At this point, it is a good idea to save the standard sweep width FID.**

LC Next

LC Thres

LC HOLD, then release

axis='p' <>

pll page <>

access the next display menu

select threshold

move threshold line to include the peak(s)

set axis to ppm

print peak frequencies in Hz and ppm

The written record of the referenced chemical shifts is necessary if the sweep width and/or transmitter offset need to be changed. If the sweep width needs to be made larger or smaller, the spectrometer will not keep a record of the reference position and the peak references will be wrong in the next spectrum. However, using the printout, the peaks can be easily re-referenced.

**Example Two:** The standard sample and the research sample have lock solvents with very different chemical shifts. The user will not be able to lock on the reference and should treat this like example three below.

**Example Three:** The standard sample has no lock solvent.  
<sup>31</sup>P research sample in CDCl<sub>3</sub> and 85% H<sub>3</sub>PO<sub>4</sub> standard.

insert the research sample

*use the standard setup for P31*

jexp1 <>	join experiment 1
LC Main Menu	activate main menu display
LC Setup	activate setup menu
LC Nucleus,Solvent	select nucleus, solvent menu
LC P31	select nucleus
LC CDCl3	select solvent
rts('current shim library') <>	retrieve shims
su <>	perform experiment set up

when set up is complete:

load='n' <>	set load shim values to no
dg <>	display acquisition parameters

lock the sample and crudely shim Z1 and Z2 and disconnect the VNMR Acquisition window  
 set nt to some reasonable number that will allow the peak(s) to be visible

ga <>	start acquisition (will wft when complete)
-------	--

When acquisition is complete:

f full <>	display full sweep width to screen
aph <>	autophase
make sure that the peak(s) are visible at this sweep width and number of transients	
LC Connect (in the VNMR acquisition window)	<u>IMPORTANT:</u>
LC Lock	If you do not turn the lock off before
LC Lock off	ejecting the sample, your reference
eject the research sample	<u>will not be</u> correct.
insert the standard sample	
<i>do not change any parameters</i>	
nt=4 <>	set number of transients
ga <>	start acquisition (will wft when complete)

When acquisition is complete:

f full <>  
 aph <>  
 LC HOLD, then release  
 RC HOLD, then release  
 LC Expand  
 LC HOLD, then release  
 nl <>  
 rl(0p) <>  
 f full <>  
 LC HOLD, then release

LC Mark  
 cr = \_\_\_\_\_

display full sweep width to screen  
 autophase  
 move first cursor to the left of reference line  
 move second cursor to the right of reference line  
 expand region inside cursors  
 move first cursor to center of line  
 select nearest line  
 set reference (in ppm)  
 display full sweep width to screen  
 move first cursor far right of spectrum, to the last  
 data point on the screen  
 request frequency  
 write down the frequency (which is in hertz)

**It is a good idea at this time to save the reference FID.**

LC Connect (in the VNMR acquisition window)

LC Lock

LC Lock off

eject the standard sample

insert the research sample

lock on the research sample without changing Z0

the research sample should now be shimmed as normal and the VNMR Acquisition window disconnected

*do not change any other parameters*

set nt so that the peak(s) will be visible

ga <>

IMPORTANT:

If you do not turn the lock off before  
 ejecting the sample, your reference  
will not be correct.

start acquisition (will wft when complete)

when acquisition is complete:

f full <>  
 aph <>  
 LC HOLD, then release

rl(cr value in hertz) <>  
 f full <>

display full sweep width to screen  
 autophase  
 move first cursor far right of spectrum, to the last  
 data point on the screen  
 set reference (in hertz)  
 display full sweep width to screen

**At this point, it is a good idea to save the standard sweep width FID.**

LC Next	access the next display menu
LC Thres	select threshold
LC HOLD, then release	move threshold line to include the peak(s)
axis='p' <>	set axis to ppm
pll page <>	print peak frequencies in ppm

The written record of the referenced chemical shifts is necessary if the sweep width and/or transmitter offset need to be changed. If the sweep width needs to be made larger or smaller, the spectrometer will not keep a record of the reference position and the peak references will be wrong in the next spectrum. However, using the printout, the peaks can be easily re-referenced.

**NOTE: If vt is used for the research sample, the standard should be run at the same temperature, if possible. Under no circumstances should samples in D<sub>2</sub>O or 85% H<sub>3</sub>PO<sub>4</sub> be run below 0°C. If the reference sample can not be run at the same temperature as the research sample, note the conditions of the reference sample and the research sample and report these in appropriate experimental sections.**